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Instructions for Use

Product name:Medical Protective Mask

Product Model: C004

Primary structure: 5 layers of protection

Executive standard: GB-19083 2010

Materials: The product consists of mask, nose clip and mask belt. The mask body has a five-layer structure composed of an outer layer, an inner layer, two filter layers and another inner layer. The inner and outer layers of the mask body are polypropylene spunbonded nonwovens, the inner layer is hot-air cotton, and the filter layer is polypropylene meltblown nonwovens; the nose clip is made of aluminum strip material; the mask is made of nylon and spandex elastic material. Products are provided in both sterile and non-sterile forms.

Executive standard: GB 19083-2010

Usage range: It is used for filtering particles in the air, blocking droplets, blood, body fluids, and secretions in the medical working environment.

Use instructions:

Step 1. Adjust and fasten the headbands to keep mask body closely to

your nose and face.

Step 2. Adjust the nose clip and press it, and keep it closely to nose and

face.

Step 3. Close to the face, and stretch the mask until it fully covers the nose, mouth and chin. Exhale sharply to check the airtightness. If air leaks around the nose, readjust the headbands and the nose clip as described in step 1 and step 2 to eliminate leakage. If air leaks at the mask edges, adjust the mask position slowly. Repeat this procedure until

the mask fits properly.

Step 4. If damaged, smelly or damp, please replace it with a new one immediately.

Warnings:

1. Before use, please ensure there is no damage or pollution on the overall appearance, good elasticity of headbands, and no damage or fracture on the nose clip. Confirm the external packaging mark, production date, and expiry date, and use it within the expiry date.
2. This product is single use only and it is strictly prohibited to reuse it. After use, it should be disposed of according to the requirements of the hospital and the environmental protection department.
3. During use, if the respiratory resistance is significantly increased (when breathing is difficult), the mask is damaged, or the mask is contaminated (when it is stained with foreign objects such as blood stains and droplets), it should be replaced in time.
4. Do not modify, clean, abuse or misuse this mask.
5. Avoid hand contact with the inside of the mask.

6. The mask shall be replaced in time when it is wet.
7. Don't use equipment in explosive atmosphere.
8. This mask does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
9. Please pay attention to avoid high temperature and humidity during storage, and keep in a clean place.

Storage conditions: The packaged products should be stored in a cool, dry, well-ventilated clean room with a temperature of $-20\text{ }^{\circ}\text{C}$ - $+30\text{ }^{\circ}\text{C}$ and a relative humidity of no more than 80%. The room should be protected from direct sunlight and away from fire, toxic and corrosive gases.

Transport conditions: The product shall be loaded and transported in covered carriage and cabin, and kept clean to avoid heavy pressure, sunshine and rain.

Manufacturer: Qingdao Hainuo Biological Engineering Co., LTD.

Address: No. 1 Haishi Hainuo New Industrial Area, Jiangshan Town Industrial Park, Laixi, Qingdao, China.

Contact information: Tel: 0532-86999302

Manufacturing Date: See package

Period of validity: 2 years