

QINGDAO HAINUO GROUP

About us

Hainuo Group was founded in 2004 and now has more than 20 subsidiary corporations. Our headquarter is located in Jiangshan industrial park in Laixi, Qingdao, which is very close to the main highway, Huashan International Golf Course, Qingdao International Airport and Qingdao Harbor. Our factory covers an area of 220,000 square meters with super geographical position and convenient transportation.

We have a dozen up-to-date factories, multiple Class 100,000 cleanrooms and Class 300,000 cleanrooms, modern logistics and storage facilities, advanced production facilities and detecting instruments. More than 2,000 professional staffs devote themselves to their job and the development of Hainuo, a lot of them are experts, professors and licensed pharmacists from Chinese Pharmaceutical Association as well as the State Council. We also have technology support and cooperation of many research institutions and companies from USA, Germany and Taiwan, China.

Hainuo keeps moving forward and we're always open to your ideas and advices. Please don't hesitate to let us know if you're interested in working with us and a deeper level of cooperation!

OUR BRANDS



HYNAUT for Overseas Market



HAINUO for Domestic Market



3强实企
STRONG FIRM
品牌质量加持

Contents Page

1. Wound Care Series

1.1 Adhesive Bandage (Base Fabric Series)	Page 3
1.2 Adhesive Bandage (Transparent Series)	Page 4
1.3 Adhesive Bandage (Childlike Series)	Page 5
1.4 Chitin Adhesive Bandage	Page 6
1.5 Alginate Adhesive Bandage	Page 7
1.6 Hydrocolloid Adhesive Bandage	Page 8
1.7 Liquid Bandage	Page 9

2. Medical Care Series

2.1 Sterile Wound Dressing	Page 10
2.2 Pressure-sensitive Adhesive Tape (Cotton Fabric Series)	Page 11
2.3 Pressure-sensitive Adhesive Tape (Non-woven & PE Series)	Page 12
2.4 Warm Patch	Page 13
2.5 Cool Patch & Fever Relief Patch	Page 14
2.6 Gauze Pad & Gauze Bandage	Page 15
2.7 Medical Disposable Mask	Page 16
2.8 Medical Nursing Pads & Paper Diapers	Page 17
2.9 Cold Pack	Page 18



3. Disinfection Series

3.1 Iodine Swabs & Cotton Swabs	Page 19
3.2 Alcohol Wipes & Cotton Balls	Page 20

4. Family Care Series

4.1 First Aid Kit	Page 21
4.2 Face Film	Page 22
4.3 Floss Pick	Page 23
4.4 Moxibustion	Page 24

5. Latex Condom Series

5. Natural Latex Condom	Page 25
-------------------------	---------

Wound Care Series

Adhesive Bandage (Base Fabric Series)

HYNAUT adhesive bandage adopts high-quality adhesive from German Henkel to ensure a super durable adhesion. Its thickened island-shaped pad has outstanding absorbency. Some of these products are water resistant and also provide a lot of flexibility and comfort. They can be used for superficial trauma or bruise in/above the shallow-layer skin, they have excellent hemostatic function and provide an ideal healing condition for wound surface and skin injury.

Features:

- ✓ German Henkel high-quality adhesive, durable adhesion
- ✓ Thickened island-shaped pad with outstanding absorbency
- ✓ Flexible and comfortable

Variety	Specification	Packing
Water Resistant & Elastic	72 x 22mm	100 pcs/box, 100 boxes/carton
Water Resistant	72 x 22mm	6 pcs/bag, 20 bags/box, 60 boxes/carton
Water Resistant	70 x 22mm	10 pcs/bag, 20 bags/box, 48 boxes/carton
Water Resistant	70 x 22mm	6 pcs/bag, 20 bags/box, 60 boxes/carton
Lightweight	70 x 18mm	12 pcs/box, 400 boxes/carton
Lightweight	70 x 18mm	26 pcs/box, 400 boxes/carton
Elastic	70 x 18mm	100 pcs/box, 100 boxes/carton
Economical	70 x 18mm	100 pcs/box, 100 boxes/carton
Ultra Large	90 x 65mm	4 pcs/box, 400 boxes/carton
Ultra Large	90 x 65mm	2 pcs/box, 400 boxes/carton
Ultra Large	70 x 60mm	4 pcs/box, 200 boxes/carton
Ultra Large	100 x 100mm	4 pcs/box, 200 boxes/carton





Your Professional Partner
of Wound Care



Adhesive Bandage (Transparent Series)

HYNAUT adhesive bandage adopts high-quality adhesive from German Henkel to ensure a super durable adhesion. The product can be used for superficial trauma or bruise in/above the shallow-layer skin, it has excellent hemostatic function and provides an ideal healing condition for wound surface and skin injury.

Features:

- ✓ Waterproof, transparent and almost invisible
- ✓ German Henkel high-quality adhesive, durable adhesion
- ✓ High elastic PU/PE material with ultrathinness
- ✓ Excellent hemostatic function



Variety	Specification	Packing
Water Resistant & Tattooed	40 x 40mm	8 pcs/box, 400 boxes/carton
Water Resistant & Tattooed	28 x 21mm	8 pcs/box, 400 boxes/carton
Elastic PE	72 x 22mm	6 pcs/bag, 40 bags/box, 30 boxes/carton
Waterproof	72 x 22mm	5 pcs/bag, 20 bags/box, 60 boxes/carton
Waterproof	72 x 22mm	6 pcs/box, 400 boxes/carton
Waterproof	62 x 25mm	10 pcs/box, 400 boxes/carton
Waterproof	69 x 28mm	10 pcs/box, 400 boxes/carton

Adhesive Bandage (Childlike Series)

HYNAUT childlike adhesive bandages are designed adorably. They all adopt high-quality adhesive from German Henkel to ensure a super durable adhesion. The products can be used for superficial trauma or bruise in/above the shallow-layer skin, it has excellent hemostatic function and provides an ideal healing condition for wound surface and skin injury.

Features:

- ✓ Adorable design, recovering childhood
- ✓ High flexibility and elasticity
- ✓ German Henkel high-quality adhesive, durable adhesion
- ✓ Excellent hemostatic function

Variety	Specification	Packing
Various Shapes Combined	90 x 65mm 72 x 38mm 60 x 45mm 76 x 38mm 72 x 22mm	18 pcs/box, 400 boxes/carton
Cartoon	70 x 18mm 56 x 19mm	5 pcs/bag, 20 bags/box, 60 boxes/carton
Cartoon	72 x 22mm	6 pcs/bag, 40 bags/box, 20 boxes/carton
Cartoon	56 x 19mm 25 x 25mm 70 x 18mm	8 pcs/box, 400 boxes/carton
Cartoon	72 x 22mm	12 pcs/box, 400 boxes/carton



Chitin Adhesive Bandage

Chitin is the only positively charged cellulin that has ever been discovered in nature, it has powerful adsorption and inhibitory effect against negatively charged harmful substances, such as bacteria, fungus and virus. Therefore, it has significant effect against staphylococcus aureus, colibacillus, Candida albicans, gram-positive bacteria etc. With histocyte-friendly glucosamine and fibrocyte suppressing, chitin can also accelerate tissue repair and limit scarring.

HYNAUT chitin adhesive bandage can be used for superficial trauma or bruise in/above the shallow-layer skin. Meanwhile it provides an ideal healing condition for wound surface and skin injury.

Features:

- ✓ Accelerates wound healing
- ✓ Inhibits scar formation
- ✓ Superior hemostatic function and pain relief
- ✓ Significant bacteria inhibition

Variety	Specification	Packing
Chitin Adhesive Bandage	72 x 22mm	8 pcs/box, 200 boxes/carton





Your Professional Partner
of Wound Care



Alginate Adhesive Bandage

Alginate is a non-toxic natural macromolecular materials, it can absorb up to 11 times its own weight in liquid. When alginate contacts with wound exudate, the calcium ion inside the absorbent pad will exchange with the sodium ion in the exudate, then forms hydrosoluble sodium alginate. It forces the exudate to enter the absorbent pad and then forms a layer of hydrogel membrane on the surface of wound. With its oxygen permeability the hydrogel membrane can prevent bacteria from entering wound, meanwhile it can promote creation of fresh skin tissue and accelerate wound healing through its moist substrate.

HYNAUT alginate adhesive bandage can be used for superficial trauma or bruise in/above the shallow-layer skin, it provides an ideal healing condition for wound surface and skin injury. This product is extremely suitable for large exudation wounds.

Features:

- ✓ Accelerates wound healing and inhibits scar formation
- ✓ Rapid hemostasis and pain relief
- ✓ High oxygen permeability and anti-bacteria
- ✓ Avoids secondary injury of wound while peeling off

Variety	Specification	Packing
Alginate Adhesive Bandage	72 x 22mm	4 pcs/box, 400 boxes/carton





Your Professional Partner
of Wound Care



Hydrocolloid Adhesive Bandage

Hydrocolloid adhesive bandage has significant effect against chronic and refractory wounds, it's also very suitable for tissue regeneration, which is very complicated and time-consuming. It provides an ideal healing condition for wound healing and is easy to remove without damaging new tissue.

HYNAUT hydrocolloid adhesive bandage can be mainly used for nursing of wounds without exudate or with slight exudate, for example superficial burn, skin abrasion, etc. Due to good sealing ability it's sterile and water resistant, even when doing sports or playing, this feature makes this product an ideal option for children.

Features:

- ✓ Larger than 3 times of water absorbency
- ✓ Sterile, water resistant, transparent, light and comfortable
- ✓ Excellent pain relief in burn treatment
- ✓ Accelerates wound healing
- ✓ No pain while peeling off, perfect for children

Variety	Specification	Packing
Hydrocolloid Adhesive Bandage	72 x 45mm	2 pcs/box, 300 boxes/carton
Hydrocolloid Adhesive Bandage	72 x 19mm	5 pcs/box, 300 boxes/carton





Liquid Bandage

Our liquid bandage doesn't include any components with pharmacological effects, it's unable to be absorbed by human bodies as well.

The liquid will quickly form a viscous and water resistant protective film on the wound surface in 10 ~ 30 seconds. This film can last about 2 days and it will play a physical barrier function.

This product can be used for nursing small wounds, bruise, incised wounds, superficial traumas and surrounding skin, etc.

Features:

- ✓ 10s~30s of rapid film formation
- ✓ 2 days of film duration
- ✓ Able to contact directly with water
- ✓ Safe and convenient to use

Variety	Specification	Packing
Liquid Bandage	10ml	1 bottle/box



Your Professional Partner
of Wound Care



Sterile Wound Dressing

HYNAUT sterile wound dressing consists of medical tape, absorbent pad and release paper. The surface is made of non-woven fabric to ensure comfort. We also adopt high-quality adhesive from German Henkel to ensure it's hypoallergenic and durable sticky. This ethylene oxide sterilized product can be used as dressing for trauma after debridement and postoperative wound, also it can be used for fixation of IV catheter.

Features:

- ✓ German Henkel high-quality adhesive, durable adhesion
- ✓ Medical grade adhesive, hypoallergenic
- ✓ Non-woven/PU surface, very comfortable
- ✓ Sterilized by ethylene oxide, safe and secure

Variety	Specification	Packing
Sterile Wound Dressing	6 x 7cm	1 pc/bag, 80 bags/box, 20 boxes/carton
Sterile Wound Dressing	7 x 9cm	1 pc/bag, 80 bags/box, 20 boxes/carton
Sterile Wound Dressing	10 x 10cm	1 pc/bag, 80 bags/box, 20 boxes/carton
Sterile Wound Dressing	10 x 15cm	1 pc/bag, 40 bags/box, 20 boxes/carton
Sterile Wound Dressing	10 x 20 cm	1 pc/bag, 40 bags/box, 20 boxes/carton
Sterile Wound Dressing	10 x 25cm	1 pc/bag, 40 bags/box, 20 boxes/carton
Sterile Wound Dressing	10 x 30cm	1 pc/bag, 40 bags/box, 20 boxes/carton
For Detaining Needle	6 x 7cm	1 pc/bag, 100 bags/box, 20 boxes/carton
I.V. Adhesive Dressing	7 x 3.5cm	100 pcs/box, 100 boxes/carton



Pressure-sensitive Adhesive Tape (Cotton Fabric Series)

HYNAUT pressure-sensitive adhesive tape adopts high-quality adhesive from German Henkel to ensure it is durable sticky and hypoallergenic. The surface is made of non-woven fabric, which is very light and easy to tear off. This product can provide durable adhesion for wound dressing, bandage, etc. so that a reliable fixation can be effected. It's used to contact directly with intact skin, instead of wound surface.

Features:

- ✓ Premium cotton fabric
- ✓ German Henkel high-quality adhesive, durable adhesion
- ✓ Medical grade adhesive, hypoallergenic
- ✓ Easy to peel off, convenient and fast

Variety	Specification	Packing
Zinc Oxide Content	2.5cm x 4m	1 roll/box, 400 boxes/carton
Cotton Fabric	0.9cm x 3 / 4 / 5m	80 boxes/carton
Cotton Fabric	0.9cm x 6 / 7 / 8 / 9 / 10m	60 boxes/carton
Cotton Fabric	2.0cm x 3 / 4 / 5m	1 roll/box, 400 boxes/carton
Cotton Fabric	2.5cm x 3 / 4 / 5m	1 roll/box, 400 boxes/carton
Cotton Fabric	26cm x 3 / 4 / 5m	60 boxes/carton
Cotton Fabric	28cm x 3 / 4 / 5m	30 boxes/carton



Pressure-sensitive Adhesive Tape (Non-woven & PE Series)

HYNAUT pressure-sensitive adhesive tape adopts high-quality adhesive from German Henkel to ensure it is durable sticky and hypoallergenic. The surface is made of non-woven fabric, which is very light and easy to tear off. This product can provide durable adhesion for wound dressing, bandage, etc. so that a reliable fixation can be effected. It's used to contact directly with intact skin, instead of wound surface.

Features:

- ✓ German Henkel high-quality adhesive, durable adhesion
- ✓ Medical grade adhesive, hypoallergenic
- ✓ Excellent breathability (PE type)
- ✓ Easy to peel off, convenient and fast

Variety	Specification	Packing
Non-woven	1.1cm x 9.1m	24 rolls/box, 40 boxes/carton
Non-woven	1.25cm x 9.1m	24 rolls/box, 40 boxes/carton
Non-woven	2.5cm x 9.1m	12 rolls/box, 40 boxes/carton
Non-woven	5.0cm x 9.1m	6 rolls/box, 40 boxes/carton
PE, breathable	1.25cm x 3 / 4 / 5m	24 rolls/box, 80 boxes/carton
PE, breathable	1.25cm x 6 / 7 / 8 / 9 / 10m	24 rolls/box, 60 boxes/carton
PE, breathable	2.5cm x 6 / 7 / 8 / 9 / 10m	12 rolls/box, 60 boxes/carton
PE, breathable	5.0cm x 6 / 7 / 8 / 9 / 10m	6 rolls/box, 60 boxes/carton
PE, breathable	1.25cm x 3 / 4 / 5 / 6 / 7 / 8 / 9 / 10m	24 rolls/box, 40 boxes/carton



Warm Patch

HYNAUT warm patch consists of highly enriched exothermic material, non-woven bag and sticker. It can be used for lower abdomen pain during menstrual period, menstrual disorder, uterus-warming for parturients and uteruswarming after abortion. Also it's suitable for soft tissue injury, arthralgia, lumbocrurol pain, stomach cold, periarthritits of shoulder, cervical spondylosis, etc. We provide you at least 8-hours continuous warming and protect your health. The plaster will start heating up after opened and it will be applied through clothes, instead of direct contact with the skin.

Features:

- ✓ Rich functions
- ✓ $\leq 70\text{ }^{\circ}\text{C}$ temperature
- ✓ ≥ 8 hrs of continuous warming
- ✓ ≤ 10 min. quick heating and easy to use
- ✓ Lightweight and comfortable



Variety	Specification	Packing
Universal type	96 x 130mm	10 pcs/plastic bag, 30 bags/carton
For neck and shoulder	54x93mm + 70x72mm + 72x70mm + 93x54mm (4 parts in one piece)	3 pcs/plastic bag, 60 bags/carton
For elbow joint and knee	90 x 135mm	5 pcs/plastic bag, 60 bags/carton
For uterus warming	96 x 130mm	5 pcs/plastic bag, 60 bags/carton
For lumbar	310 x 125mm	3 pcs/plastic bag, 60 bags/carton
For feet	70 x 90mm	6 pcs/plastic bag, 60 bags/carton

Cool Patch & Fever Relief Patch

HYNAUT cold compress plaster consists of non-woven backing layer, gel layer and PE membrane layer. It is used as physical cooling and cold therapy for the whole body and only the closed soft tissue. You can also tailor the plaster into appropriate size according to different injured parts. Every piece of plaster can continuously take effect for 4~8 hours.

Features:

- ✓ Quick cooling and outstanding cooling effect
- ✓ Easy to carry and use
- ✓ Suitable for all closed soft tissues

Variety	Specification	Packing
Cool Patch	80 x 120mm	2 pcs/bag, 2 bags/box, 300 boxes/carton
Fever Relief Patch (Adult)	50 x 120mm	2 pcs/box, 300 boxes/carton
Fever Relief Patch (Child)	40 x 110mm	3 pcs/box, 300 boxes/carton
Fever Relief Patch (Apple Smell)	40 x 110mm	4 pcs/box, 300 boxes/carton
Fever Relief Patch (Strawberry Smell)	40 x 110mm	4 pcs/box, 300 boxes/carton
Fever Relief Patch (Orange Smell)	40 x 110mm	4 pcs/box, 300 boxes/carton
Fever Relief Patch (Mixed Fruit)	40 x 110mm	4 pcs/box, 300 boxes/carton



Gauze Pad & Gauze Bandage

HYNAUT gauze pad and gauze bandage are made of 100% absorbent cotton, some models also include barium line, which is detectable by X-ray. It is sterilized by ethylene oxide and can be used for protection and cleaning of skin, mucosa and wound, during surgery it can be used for absorption of exudate as well. These products are very convenient essentials for hospital, community health service center and clinic.

Features:

- ✓ 100% cotton, excellent absorbency
- ✓ 8 gauze layers, great protection for wound
- ✓ Sterilized by ethylene oxide, safe and secure
- ✓ Soft and comfortable



Variety	Specification	Packing
Sterile Gauze Pad	7.5 x 7.5cm	2 pcs/bag, 50 bags/box, 20 boxes/carton
Sterile Gauze Pad	8 x 10cm	5 pcs/bag, 25 bags/box, 20 boxes/carton
Gauze Bandage	6cm x 6m	1 roll/bag, 600 bags/carton
Gauze Bandage	8cm x 6m	10 rolls/bag, 120 bags/carton

Medical Disposable Mask

HYNAUT disposable medical mask is an ear-hook mask with high safety through ethylene oxide sterilization.

There're 3 layers of construction: internal layer is made from PP non-woven and so is the outer layer, middle layer is made from PP melt-blown fabric. They can effectively prevent pathogenic microorganisms, body fluid and particulates from penetrating physical barrier directly.

This product can be used for covering the mouth, nose and jaw for clinical staffs during noninvasive process. It's also suitable for disposable hygiene care in common environment.



Features:

- ✓ Multi-layer prevention and protection
- ✓ Bacterial filtration efficiency $\geq 95\%$
- ✓ Sterilized by ethylene oxide, safe and secure

Variety	Specification	Packing
Ear-hook, flat type	17.5 x 9.5cm	50 pcs/bag/box, 40 boxes/carton
Ear-hook, flat type	17.5 x 9.5cm	3 pcs/bag, 20 bag/box, 40 boxes/carton
Ear-hook, flat type	17.5 x 9.5cm	5 pcs/bag, 20 bag/box, 40 boxes/carton
Ear-hook, flat type	17.5 x 9.5cm	10 pcs/bag, 1 bag/box, 200 boxes/carton
Ear-hook, folding type	15.5 x 14.5cm	1 pc/bag, 3 bag/box, 200 boxes/carton
Ear-hook, folding type	15.5 x 14.5cm	1 pc/bag, 5 bag/box, 200 boxes/carton

Medical Nursing Pads & Paper Diapers

Variety	Specification	Packing
Nursing pads for adult	60 x 90cm	10 pcs/bag, 30 bags/carton
Nursing pads for adult	60 x 90cm	22 pcs/bag, 12 bags/carton
Nursing pads for adult	60 x 60cm	26 pcs/bag, 18 bags/carton
Paper diapers for adult	60 x 80cm	10 pcs/bag, 8 bags/carton
Nursing pads for adult	60 x 90cm	10 pcs/bag, 16 bags/carton
Nursing pads for adult	60 x 60cm	15 pcs/bag, 16 bags/carton
T-shape baby diapers	25 x 31 x 10cm	24 pcs/bag, 24 bags/carton
Nursing pads for baby	25 x 35cm	80 pcs/bag, 12 bags/carton
Anti-galactorrhea nursing pads	12 x 13cm	30 pcs/bag, 24 bags/carton
Nursing pads for puerpera	30 x 60cm	1000 g/bag, 20 bags/carton
Childbed (cotton)	60 x 90cm	6 pcs/bag, 24 bags/carton
Childbed (non-woven)	60 x 90cm	10 pcs/bag, 24 bags/carton
Nursing pads for puerpera (sanitary towel type)	42 x 18cm	6 pcs/bag, 30 bags/carton



Instant Cold Pack

Features:

- ✓ Activates easily to provide instant cold
- ✓ Instant aches and pain relief
- ✓ Provides an anesthetic effect
- ✓ First aid kit




Iodine Swabs & Cotton Swabs

HYNAUT medical cotton swabs are made of medical absorbent cottons and fine bamboo sticks. We select 100% natural cotton and we provide sufficient quantity on each swab. This ethylene oxide sterilized product can be used for absorption of disinfectant and the disinfection of intact skin. It's a very convenient essential for hospital, community health service center, clinic and household.

Features:

- ✓ 100% natural cotton
- ✓ Fine bamboo stick
- ✓ Sufficient quantity on each single head
- ✓ Sterilized by ethylene oxide, safe and secure



Variety	Specification	Packing
<i>Iodine Swabs</i>	<i>0.15ml</i>	<i>24 pcs/box, 600 boxes/carton</i>
<i>Iodine Swabs (paper box)</i>	<i>0.15ml</i>	<i>20 pcs/box, 400 boxes/carton</i>
<i>Iodine Swabs (PVC box)</i>	<i>0.15ml</i>	<i>20 pcs/box, 400 boxes/carton</i>
<i>Sterile, single head</i>	<i>8cm</i>	<i>100 pcs/bag, 600 bags/carton</i>
<i>Single head</i>	<i>10cm</i>	<i>50 pcs/bag, 960 bags/carton</i>
<i>Double heads, bamboo stick</i>	<i>8cm</i>	<i>100 pcs/bag, 480 bags/carton</i>
<i>Double heads, PVC stick</i>	<i>8cm</i>	<i>100 pcs/bag, 192 bags/carton</i>
<i>Double heads, paper stick</i>	<i>8cm</i>	<i>100 pcs/bag, 192 bags/carton</i>



Alcohol Wipes & Cotton Balls

Variety	Specification	Packing
Alcohol Wipes	6 x 3cm	50 pcs/box, 200 boxes/carton
Alcohol Wipes	6 x 6cm	100 pcs/box, 100 boxes/carton
Alcohol Wipes	12 x 12cm	50 pcs/box, 60 boxes/carton
Cotton Balls	10g	400 bags/carton
Cotton Balls	50g	100 bags/carton
Iodine Cotton Balls	25 grains/bottle	240 bottles/carton
Iodine Cotton Balls	30 grains/bottle	240 bottles/carton
Iodine Cotton Balls	30 grains/bottle	240 bottles/carton
Alcohol Cotton Balls	25 grains/bottle	240 bottles/carton
Alcohol Cotton Balls	30 grains/bottle	240 bottles/carton
Alcohol Cotton Balls	30 grains/bottle	240 bottles/carton
Alcohol Cotton Balls	50 grains/bottle	180 bottles/carton
Alcohol Cotton Balls	90 grains/bottle	108 bottles/carton

Family Care Series


 Your Professional Partner
of Wound Care





First Aid Kit



Hynaut provides you various options for customizing your own first aid kits, for example home first aid kit, outdoor first aid kit, first aid kit for vehicles, earthquake first aid kit, etc. You can choose the products from both the product product portfolio of our own and from other suppliers, then put the products into the kit and make it a whole new commodity.

Variety	Packing
Camo	30 bags/carton
Red	60 bags/carton
Camo Waist Bag	40 bags/carton



Your Professional Partner
of Wound Care




Face Film

Hynaut face film can help skin to replenish moisture, it whiten skin, inhibit melanin production, brighten facial complexion and dilute spots. It will help to exclude the skin metabolism of the skin and the accumulation of oils and substances. The moisture in the mask penetrates into the stratum corneum of the epidermis, the skin becomes soft, and the skin is naturally bright and elastic.



Variety	Specification	Packing
Face Film	25ml/pc	5 pcs/box, 44 boxes/carton



HCG Test Strip

HYNAUT HCG test strip uses colloidal gold method, which is a simple and reliable method for testing early pregnancy. It doesn't require any special equipment and reagent, and it can provide an intuitional and precise diagnosis in a very short period of time.

Variety	Packing
Pen Type	1 pc/box, 480 boxes/carton
Bar Type	100 pcs/box, 50 boxes/carton
Card Type	1 pc/box, 480 boxes/carton

Floss Pick

HYNAUT floss picks are very strong and durable, it has no fraying during use.

Features:

- ✓ Strong and durable
- ✓ Mint fragrance
- ✓ Multifunctional

Variety	Packing
Floss picks	12 pcs/box, 1000 boxes/carton
Floss picks	20 pcs/box, 400 boxes/carton
Floss picks	60 pcs/box, 200 boxes/carton
Floss picks	50 pcs/bag, 300 bags/carton
Floss picks	32 pcs/bag, 288 bags/carton
Floss picks	50 pcs/bag, 384 bags/carton
Floss picks	50 pcs/box, 288 boxes/carton





Moxibustion

Moxibustion is also called moxibustion therapy, which is an very important part of Acupuncture & Moxibustion in Chinese Medicine. The use of moxibustion therapy originates from the stone age, in which the human beings just mastered fire skills. The Chinese ancestors made moxa leaf into moxibustion materials, then used the generated heat to stimulate acupuncture points or specific parts in body surface. The disordered function of physiology and biochemistry can be adjusted through stimulating the activity of Channel-qi, and the disease preventing and curing can be realized as well.

Effects & Features:

- ✓ Soothes the nerves, helps sleep, relieves pressure
- ✓ Expels wind-damp, dissipates cold, strengthens body resistance and eliminates etiological factors
- ✓ Warms and activates meridians
- ✓ Relieves palace cold and menstrual cramps
- ✓ Nurses spleen and stomach, improves digestion
- ✓ Protection against the cold, boosts immunity

Variety	Specification	Packing
Smokeless Moxa-Stick	Smokeless	8 pcs/box, 100 boxes/carton
Pure Moxa-Stick	3 / 5 / 8 / 10 Years Aged Moxa	10 pcs/box, 100 boxes/carton
Pure Moxa-Cone	For Women / 5 Years Aged Moxa	54 pcs/box, 100 boxes/carton
Wormwood Packet for Foot Bath	10g/pack	20 packs/box, 60 boxes/carton
Moxibustion Plaster	Diameter 7mm, Height: 40mm	30 pcs/box, 30 boxes/carton




Latex Condom Series



Your Professional Partner
of Wound Care



Natural Latex Condom



Variety	Width	Length	Thickness	Style	Capacity of Lubricant	Type of Lubricant	Odour	Packaging
Graphene Reinforced Condom	52mm ± 2mm	160mm / 180mm	0.03mm / 0.04mm / 0.05~0.06mm	Glossy / Ribbed / Dotted / Ribbed + Dotted	300 / 500 / 700 / 1000 mg	Silicone Oil / Water-Based	Fruity / Potpourri / Odorlessness	Aluminum Foil Bag / Butter Box
Latex Condom	52mm ± 2mm	160mm / 180mm	0.03mm / 0.04mm / 0.05~0.06mm	Glossy / Ribbed / Dotted / Ribbed + Dotted	300 / 500 / 700 / 1000 mg	Silicone Oil / Water-Based	Fruity / Potpourri / Odorlessness	Aluminum Foil Bag / Butter Box

Contact us!



QINGDAO HAINUO GROUP

Headquarter Address:

*Hainuo Building, No.1 Guangdong Road,
Jiangshan International Industrial Park,
Laixi City, Qingdao City, Shandong
Province, China*

International Department Address:

*Room 501, Building 4, No.187 Jinshui
Road, Licang District, Qingdao City,
Shandong Province, China*

Tel.:

+86-532-84670375

Website:

www.hynaut.com.cn

www.hainuocn.com

E-mail:

wangcheng@hynaut.cn