

HYNAUT
海氏海诺®

Hynaut Disposal Vinyl Gloves



OUR BRANDS



HYNAUT for Overseas Market



HAINUO for Domestic Market

Qingdao Hainuo Biological Engineering Co., Ltd.
Add: NO.1 Guangdong Road, Jiangshan Industrial Zone Laixi City, Qingdao, China
Tel: +86 532 86999302 Website: <https://www.hynautgroup.cn/>

Applications

- ✓ Medical Examinations in Assisted Living Facilities
- ✓ Medical Examinations & Non-Surgical Procedures in Medical Facilities
- ✓ Home Health Care Services
- ✓ Day cares

Features

- ✓ Low cost.
- ✓ Smooth inner and outer surface.
- ✓ High tactile sensitivity.
- ✓ Not made with natural rubber latex.
- ✓ Ambidextrous fitting.

Size Available

Small(S)	Hand width 70-80 mm , 3.1 ± 0.1 grams
Medium(M)	Hand width 80-90 mm , 3.5 ± 0.1 grams
Large (L)	Hand width 90-100mm , 3.9 ± 0.1 grams
X Large (XL)	Hand width 100-110 mm , 4.3 ± 0.1 grams

Color

Clear	
-------	---

Packing

100pcs/box 10boxes/carton

Box	240 × 125 × 63mm
-----	------------------

Loading Capabilities

20 GP	1500 CTNs
40 GP	3100 CTNs
40 HP	3450 CTNs

Approvals

EN455

Cautions

- Single use only.

Qingdao Hainuo Biological Engineering Co., Ltd.

Add: NO.1 Guangdong Road, Jiangshan Industrial Zone Laixi City, Qingdao, China

Tel: +86 532 86999302 Website: <https://www.hynautgroup.cn/>

- Avoid direct sunlight.

Storage

Store in dry and ventilated environment, room temperature should not be above 37°C.

Shelf life

5 years.