



Hynaut Alcohol Wipes

FDA Certification

Establishment Registration & Device Listing

[FDA Home](#) [Medical Devices](#) [Databases](#)



1 result found for Establishment Registration
or FEI Number : 3008808560

[New Search](#)

Establishment Name	Registration Number	Current Registration Yr
QINGDAO HAINUO BIOLOGICAL ENGINEERING CO., LTD. CHINA	3008808560	2020
Floss, Dental - Floss Pick		Manufacturer
First Aid Kit Without Drug - First Aid Kit		Manufacturer
Pad, Alcohol, Device Disinfectant - Alcohol Wipes, Alcohol Prep Pad		Manufacturer
Dressing, Wound, Hydrophilic - Wound Dressing, I.V. Adhesive Dressing		Manufacturer
Tape And Bandage, Adhesive - Medical Tape, Adhesive Bandage, Wound Plaster		Manufacturer
Face Mask (Except N95 Respirator) For General Public/Healthcare Personnel Per lie Guidance		Manufacturer
Mask, Scavenging - Disposable Mask		Manufacturer
Disinfectant, Medical Devices - Iodine Swab		Manufacturer
Pack, Hot Or Cold, Disposable - Warm Patch, Cold Pack, Cool Patch, Fever Relief Patch, Cool Gel, Freeze Gel		Manufacturer
Bandage, Elastic - Elastic Bandage, Gauze Bandage		Manufacturer
Non-Surgical Isolation Gown		Manufacturer
Gauze / Sponge, Nonresorbable For External Use - Gauze Pad, Underpad, Cotton Balls		Manufacturer
Cover, Shoe, Operating-Room		Manufacturer
Mouthguard, Over-The-Counter - Face Mask		Manufacturer
Applicator, Absorbent Tipped, Non-Sterile - Cotton Swab		Manufacturer
Cap, Surgical		Manufacturer

MSDS

FJT

Page 1 / 10
MSDS NO.: A001R2005190402
Creation Date: May 22, 2020

Applicant: QINGDAO HAINUO BIOLOGICAL ENGINEERING CO., LTD.
Address: NO.1 GUANGDONG ROAD, JIANGSHAN TOWN, LAIXI, QINGDAO, CHINA

Sample information:

Sample name: Alcohol Wipes
Model: /
Sample composition/raw material (supplied by client): Refer to section 3 of the MSDS"Composition/Information on Ingredient"
Edit period: From May 19, 2020 to May 22, 2020

Required service : As specified by the client, to edit Material Safety Data Sheet (MSDS) by the submitted sample information.

Summary : As specified by the client, This safety data sheet was prepared in accordance with UN GHS Rev.8, The EU CLP REGULATION (EC) No 1272/2008, and US OSHA Hazard Communication Standards (29 CFR 1910.1200). Please refer to the attached report for details.

Shanghai Fajin Testing Technology Co.,Ltd



Precautionary statement(s)

Prevention precautionary statements: P210: Keep away from heat/sparks/open flames/hot surfaces – No smoking

Response Precautionary Statements: P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing
P337+313: If eye irritation persists get medical advice/attention

Storage Precautionary Statements: P403+235: Store in a well ventilated place. Keep cool

Disposal Precautionary Statements: P501: Dispose of contents and container in accordance with local / regional / national / international regulations.

2.3 Description of any hazards not otherwise classified

None. See section 11 for more detailed information on health effects and symptoms.

Section 3 - Composition/Information on Ingredient

Substance () Preparation (√) Article ()

Composition: Non-woven cloth saturated with the following liquid.

Chemical name	CAS No.	Content (V/V)
75% Ethanol	64-17-5	50%
RO Pure water	7732-18-5	48.5%
Double-chain quaternary ammonium salt	-	0.15%

Abbreviation: CAS No. is Chemical Abstract Service Registry Number.

Section 4 - First Aid Measures

4.1 Description of first aid measures:

Eye Contact: Flush opened eyes with plenty of water for several minutes. If eye irritation persists, get medical advice/attention

Skin Contact: No need for first aid is anticipated. Wash off with water if necessary.

Inhalation: No need for first aid is anticipated under normal usage (single exposure of small amount) .

Ingestion: Accidental ingestion at a level high enough to be dangerous to health is unlikely. Rinse mouth with water. Obtain medical attention if symptoms appear or if large quantities have been ingested.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labeling (see section 2.2) and/or in section 11.

4.3 Indication of any immediate medical attention and special treatment needed

Directions



1. Open The Cover



2. Pull Out A Wipe



3. Pops Up The Wipe



**4. Close Lid To
Retail Moisture**

Model	Wipe Size	Barrel Size	Weight	Carton Size
60pcs	10*15cm	8*12.5cm	160g	50*33*26cm
70pcs	15*18cm	10*12.5cm	180g	53*42.5*27cm
100pcs	10*15cm	10*17cm	260g	53*42.5*27cm

青岛海诺生物工程有限公司

QINGDAO HAINUO BIOLOGICAL ENGINEERING CO., LTD.

中国山东省青岛市莱西姜山外商工业园广东路1号海诺大厦
NO.1 GUANGDONG ROAD, JIANGSHAN TOWN, LAIXI,
QINGDAO, CHINA 266603
T: 0532-86999302
E: info@hynaut.cn
U: www.hynautgroup.cn